

# VAPCO LINIMENT

**For Equines and Dogs**

**For External use only**

## COMPOSITION: Each litre contains:

|         |       |
|---------|-------|
| Menthol | 16 gm |
| Camphor | 33 gm |

## PROPERTIES:

The liniment acts through its components, as a rubefacient and counter-irritant to the skin, stimulating and increasing the local circulation; and due to the dilatation of the local blood vessels, a sensation of coldness is produced followed by analgesic effect and relief of itching.

## INDICATIONS:

- Alleviation of swellings and treatment of myositis, bursitis, sprain tendons, tendovaginitis, tendinitis, lameness, bruises, colds and rheumatism.
- Relief of muscular soreness, aches and stiffness; and arthritic pains due to over-exertion, fatigue, strenuous exercise and competitive activity.

## APPLICATION:

- **Lameness, sprains, swellings etc...:** Foment the affected area with hot water, then rub well with the liniment 2-3 times daily.
- **Rheumatism and colds:** Foment the affected area (chest, neck, shoulder, back etc...) with hot water, then rub well with the “**Spirit lotion**” ; then cover, at night, with heavy blanket and protect the animal from drafts.
- **For blistering:** Rub the targeted part with liniment, then cover with gauze saturated with the “**Spirit lotion**”.  
Repeat twice daily, if necessary.
  - **Preparation of the spirit lotion:**  
Mix(60 ml liniment + 60 ml water + 300 ml methylated spirit).
- **Setting up and leg brace in equines for racing:** Mix liniment, water, vinegar (1:10:4 parts) respectively, then rub the legs thoroughly and bandage untightly.
- **For cooling body wash in equines:** Mix liniment, water, vinegar (1:6:2 parts) respectively, to be used as a therapeutic body wash when refreshing and coolant is indicated for tired, over-heated horse.

## PRECAUTIONS:

- Do not apply to irritated or broken skin; so abrasions and wounds should be treated until healing before applying liniment.
- Avoid contact with eyes.

## STORAGE:

Keep in a dry, dark place below 30 °C.

**PACKING:** (100, 250, 500, 1000) ML.